

Muffle Furnace

The testronix muffle furnace offers high accuracy and has application in varied industries.

The instrument is used to calculate the ash content of the samples. The properties of any substance can be analyzed on the basis of percentage of ash content the material contains.



Features

- Highly accurate test results under uniform temperatures.
- Material Mild Steel
- Insulated door
- High Grade density imported Glass Wool Blanket insulation.
- Maximum Thermal Efficiency
- Digital timer for better data readability
- Accurate results
- User friendly and easy to operate

Key specifications

- Inner Chamber Size- • 4 x 4 x 9 Inch
- Display 7 Segment LED display
- Accuracy $\pm 5^{\circ}\text{C}$
- Least Count/Resolution 1°C
- Power 20A,220VAC
- Timer Range Upto 999 hrs
- Temperature Range Ambient to 900°C , 1200 C & 1400 C
- Insulation High grade ceramic insulation

