

Specific Gravity Balance Digital – Advanced model

The specific gravity balance is one of the most important instruments. It provides a rough assessment of the physical properties. The device provides the specific gravity on the digital screen accurately. The instruments provides accurate and precise reading and is easy to operate.



Product specifications:

- Max weight capacity 200gm
- Least count 0.01gm
- Pan size 120mm
- Specific gravity range 1.0 to 3.0
- Calibration 200gm brass weight

Product features:

- LCD display for accurate readings
- AC adapter power supply
- Steady Hanging hook and beaker platform

