

## Specific Gravity Balance Digital – Advanced model

The specific gravity balance is one of the most important instruments. It provides a rough assessment of the physical properties. The device provides the specific gravity on the digital screen accurately. The instruments provides accurate and precise reading and is easy to operate.



### Features

- LCD display for accurate readings
- AC adapter power supply
- Steady Hanging hook and beaker platform

### Specifications

- Max weight capacity 200gm
- Least count 0.01gm
- Pan size 120mm
- Specific gravity range 1.0 to 3.0
- Calibration 200gm brass weight

