



**TESTRONIX<sup>®</sup>**  
INSTRUMENTS




## Base Clearance Unit Membrane

This is a precise instrument which helps in measuring the depth of bottles. Base clearance membrane and value clearing unit is a high accurate instrument which consists of a digital gauge to provide the user with accurate test results. It is a simple to use and easy to operate instrument. The digital micron gauge provides consistent test results

It complies with international test standard and is a handy instrument for bottle depth testing.

**MODEL NO. TX- BCVC**



Conical Tip



Digital Dial Gauge





## FEATURES

- Easy and accurate readout
- Comes with Zero Tolerance.
- Precise measurement operations with flat base surface plate.
- Repeatability in test measurements of bottles base

## SPECIFICATION

Range	0 – 12.7mm
Least count	0.01mm
Accuracy	±0.1 mm
Display	Digital

Thank you customers for choosing us as your partners in growth !



\*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

## OUR OTHER PRODUCTS :-

Vacuum Leak Tester Touch Screen	Torque Tester Touch Screen	Hot Water Bath	Bottle Burst Tester Digital	Polariscope	Flip Cap opening tester Touch Screen

