



Laminar Air Flow - LFI 100

Testronix laminar airflow is a highly useful instrument utilized in testing laboratories to create and maintain a clean working environment. It serves as a protective barrier against external contaminants that could potentially harm the samples or compromise the accuracy of test results.

One notable feature of the Testronix laminar airflow is its side panels, which are constructed using transparent Plexiglass frames. This design allows for easy observation and monitoring of the work area without compromising the integrity of the clean environment.

To enhance visibility and facilitate the testing process, the instrument is equipped with fluorescent illuminated lights at the workstation. This ensures that the workspace is well-lit and that operators can clearly observe and manipulate samples.

This quality testing instrument incorporates a HEPA (High-Efficiency Particulate Air) filter in the work area. HEPA filters are known for their exceptional ability to trap and remove airborne particles, ensuring a high level of cleanliness and preventing contaminants from entering the workspace.

MODEL NO. TX - LFI 100





FEATURES

- Side panels of the instruments are designed using transparent Plexi glass frames.
- The instrument is duly fitted with fluorescent illuminated lights at the workstation.
- Equipped with HEPA filter in the work area
- Standardized and controlled lightening conditions
- Easy to use Menu driven system
- Wide Viewing Area
- High-Quality Light sources used
- Plexi-glass frames for sample viewing

SPECIFICATION

- Power Supply: 220 V , 50 Hz, Single Phase.
- Filter: HEPA filter
- Lights: FL Fluorescent Light
- UV Blacklight (Ultra Violet Black Light)
- Blower Motor: 1/5 HP
- Blower Timer: Hrs.: min
- Viewing Area Material: Transparent Plexiglass
- No. of UV lights: 1
- No. of FL lights: 1
- Main Body Material: Steel

Thank you customers for choosing us as your partners in growth !



*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

OUR OTHER PRODUCTS :-

