



**TESTRONIX<sup>®</sup>**  
INSTRUMENTS



## Notch Cutter

The notch cutter has immense importance for making a notch in the izod charpy samples. This is specimen preparation equipment that makes samples for izod dart impact tester to test the plastic sheets. This equipment assists in making a notch (upto 5mm). The Notch cutter is durable and reliable instrument. It is a very easy to operate and convenient instrument. The instrument consist of a motorized mechanism which further makes the operation easy.

**MODEL NO. TX-NC100**





## FEATURES

- Highly accurate Cutting for sample
- Strong gripping fixtures for Izod V-Notch sample preparation
- Cutting system motorised operation
- V-Notch measurement thought precession micrometer Top View Front View Side View Accessories
- Notch depth comparator feature
- Robust & high strength body

## SPECIFICATION

- Angle of Notch 45 degree
- Micrometer Range 0 - 25mm Micrometer
- Least Count 0.01mm
- Depth of Notch 2.5mm  $\pm$  0.05
- Power 220V Ac, 5 Amp., Single phase, 50 Hz
- Dimensions: 314 x 359 x 237 mm
- Mild steel structure

Thank you customers for choosing us as your partners in growth !



\*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

## OUR OTHER PRODUCTS :-

					
MFI-Neo 2.0	Izod Charpy Impact Tester Touch Screen	Peel Strength Tester Digital	Tearing Strength Tester Panorama	Heat Seal Tester	Hot Air Oven Digital

