



## Fume hood - FH 100

Testronix Laboratory fume hoods are specially designed to ensure the safety of laboratory personnel by effectively preventing the escape of contaminants, such as chemical vapors, dust, mists, and fumes into the laboratory environment. This lab testing instrument will provide a reliable containment solution that safeguards both the experiment and the personnel involved.

The primary purpose of Testronix Laboratory fume hoods is to create a controlled airflow system that captures and removes hazardous substances generated during laboratory procedures. It features a robust construction with a durable frame and a transparent front sash, allowing for easy visibility and accessibility to the work area.

It is equipped with high-quality fans and filters that effectively capture and filter out harmful contaminants, ensuring the purity of the laboratory air. The fans create a negative pressure within the hood, which pulls air and contaminants away from the operator, preventing their dispersion into the laboratory space.

Testronix Laboratory fume hoods also incorporate various safety features, including adjustable sash heights, ergonomic design, and proper lighting, to enhance user comfort and facilitate efficient operations. The adjustable sash height allows users to optimize the airflow and containment while maintaining a clear view of the work area. The ergonomic design ensures comfortable working conditions, reducing fatigue and improving productivity.

### MODEL NO. TX - FH 100





## FEATURES

- Ergonomic and Energy Efficient Design
- Fumes are completely eradicated from the workplace.
- Thick glass covers in front for quiet and safe operation.
- Excellent chemical resistance and high-performance airflow control.
- Low maintenance.
- Quiet operation, because the extract fan is some distance from the operator

## SPECIFICATION

- Power Supply: 230V, 50 Hz, Single Phase AC
- Motor: ¼ HP, Single Phase, 1440 Rpm
- Fan Blower: 4 inches
- Body: Stainless Steel Chamber
- Size: 2 feet x 2 feet

*Thank you customers for choosing us as your partners in growth !*



*\*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.*

## OUR OTHER PRODUCTS :-

